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## References

1. Veteran's Administration Central Office. (2006 February 8). Cessation of topical spray benzocaine usage to anesthetize the surfaces of the nasopharynx, oropharynx, laryngotracheal region and airway. affected products include, but are not limited to: • Hurracaine® spray (benzocaine 20%) • Cetacaine® spray (benzocaine 14%, butyl aminobenzoate 2% and tetracaine 2%) • Topex® spray (benzocaine 20%). Retrieved August 21, 2007 from <http://www.va.gov/ncps/alerts/Benzocaine-WWW.pdf>
2. Food and Drug Administration. (2006, February 10). Benzocaine sprays marketed under different names, including Hurracaine, Topex, and Cetacaine. Retrieved August 21, 2007 from <http://www.fda.gov/cder/drug/advisory/benzocaine.htm>
3. Astra-Zeneca Pharmaceuticals. (2005, May). EMLA Cream. Retrieved August 21 from <http://astrazeneca-us.com/pi/emla.pdf>
4. Wagner A, Marosi C, Binder M, Roggla G, Staudinger T, Keil F, Locker GJ, Frass M. Fatal poisoning due to dapsone in a patient with grossly elevated methaemoglobin levels. Br J Dermatol. 1995 Nov;133(5):816-7.
5. Yang JJ, Lin N, Lv R, Sun J, Zhao F, Zhang J, Xu JG. Methemoglobinemia misdiagnosed as ruptured ectopic pregnancy. Acta Anaesthesiol Scand. 2005 Apr;49(4):586-8.
6. National Academy of Sciences - National Research Council Academy of Life Sciences. (1981). The health effects of nitrate, nitrite and N-nitroso compounds. Washington DC: National Academy Press.
7. Taylor MB, Christian KG, Patel N, Churchwell KB. Methemoglobinemia: Toxicity of inhaled nitric oxide therapy. Pediatr Crit Care Med. 2001 Jan;2(1):99-101.
8. Gupta SK, Gupta RC, Seth AK, Gupta AB, Bassin JK, Gupta A. Methaemoglobinaemia in areas with high nitrate concentration in drinking water. Natl Med J India. 2000 Mar-Apr;13(2):58-61.

9. Murone AJ, Stucki P, Roback MG, Gehri M. Severe methemoglobinemia due to food intoxication in infants. *Pediatr Emerg Care*. 2005 Aug;21(8):536-8.
10. Kross BC, Ayebo AD, Fuortes LJ. Methemoglobinemia: nitrate toxicity in rural America. *Am Fam Physician*. 1992 Jul;46(1):183-8.
11. Dusdieker LB, Getchell JP, Liarakos TM, Hausler WJ, Dungy CI. Nitrate in baby foods. adding to the nitrate mosaic. *Arch Pediatr Adolesc Med*. 1994 May;148(5):490-4.
12. Hanukoglu A, Dan on PN. Endogenous methemoglobinemia associated with diarrheal disease in infancy. *J Pediatr Gastroenterol Nutr*. 1996 Jul;23(1):1-7.
13. Ohashi K, Yukioka H, Hayashi M, Asada A. Elevated methemoglobin in patients with sepsis. *Acta Anaesthesiol Scand*. 1998 Jul;42(6):713-6.
14. Preiser JC, Reper P, Vlasselaer D, Vray B, Zhang H, Metz G, Vanderkelen A, Vincent JL. Nitric oxide production is increased in patients after burn injury. *J Trauma*. 1996 Mar;40(3):368-71.
15. Ash-Bernal R, Wise R, Wright SM. Acquired methemoglobinemia: a retrospective series of 138 cases at 2 teaching hospitals. *Medicine*. 2004 Sep;83(5):265-73.
16. Bilgin H, Ozcan B, Bilgin T. Methemoglobinemia induced by methylene blue perturbation during laparoscopy. *Acta Anaesthesiol Scand*. 1998 May;42(5):594-5.
17. Fitzsimons MG, Gaudette RR, Hurford WE. Critical rebound methemoglobinemia after methylene blue treatment: case report. *Pharmacotherapy*. 2004 Apr;24(4):538-40.
18. Barker S, Curry J, Morgan S, Bauder W. (2006). New pulse oximeter measures methemoglobin levels in human volunteers. Paper presented at the annual meeting of the Society for Technology in Anesthesia, San Diego, CA.
19. Barker SJ, Tremper KK, Hyatt J. Effects of methemoglobinemia on pulse oximetry and mixed venous oximetry. *Anesthesiology*. 1989 Jan;70(1):112-7.

Instruments containing Masimo Rainbow SET technology are identified with the Masimo Rainbow SET logo. Look for the Masimo designation on both the sensors and monitors to ensure accurate monitoring when needed most

Masimo Corporation 40 Parker Irvine, California 92618 Tel 949-297-7000 Fax 949-297-7001 [www.masimo.com](http://www.masimo.com)

© 2007 Masimo Corporation. All rights reserved. Masimo, Rad, SET,  and Signal Extraction Technology are registered trademarks of Masimo Corporation. Pulse CO-Oximeter and Rad-57 are trademarks of Masimo Corporation. Rainbow, SpCO, CoHB, MethB are trademarks of Masimo Laboratories.